


Our autoclave features -  
ALL STANDARD!



Reduces cooling times

Cooling fan

Essential for liquid loads

Easy to use

Single action lid closure

Time saving

Protects against boil-over

Exhaust speed control

Ensures sterilisation

Protects against overcooking

Monitors load temperature, not chamber temperature

Protects against boil-over

Load sensed process timing

Uses armoured probe cable

Self-contained

Inbuilt water reservoir

Inbuilt condensate tank

No drain connections

No water connections

Additional safety

Patent 8-point lid locking

Longer chamber and lid life

Cycle start up delay

Process ends when required

Fo calculation

Clear indication of pathogen kill rate

Easy to relocate

True mobility

Mounted on castors

Advanced media programs

Media melting

Media warming

316L stainless steel chamber

Longer service life

USB interface

Easy cycle data access

Panel mounted printer

Convenient record keeping

Cycle verification

Advanced controller

60 programmes

Clear graphic display

Intuitive and easy to use

Handy user prompts throughout