

PRODUCT DATA SHEET



CV0001

Screw Neck Vial, Tall Form, With Unattached PP Closure



Product comes with unattached PP Screw Cap

Specifications

Vial Material	Glass Borosilicate 3 (Autoclavable)
Capacity	14ml
Outer Diameter	20mm
Height	72mm
Closure Type	Screw, Black Polypropylene
Closure Liner	Aluminium
Pack Quantity	100
Packaging	Cardboard Box & Plastic Bag
Pack Length	303mm
Pack Width	202mm
Pack Height	79mm
Pack Weight	1.5Kg
Case Quantity	6packs (600vials)
Case Length	440mm
Case Width	320mm
Case Height	260mm
Case Weight	11.1Kg
Pallet Quantity	204packs (20400Vials)

What's in the box?



Streamline your workflow and simplify your storage with our pack of 100 vials and caps. They neatly stack and nest together, maximizing your space and allowing convenient access whenever needed.

Unrivalled Quality

Each vial is crafted using premium-quality, chemically-resistant glass that ensures excellent durability and protection, and providing unparalleled clarity.

Range of Sizes

We understand that every project is unique. That's why our new range is available in a variety of sizes, from 2 ml to 100 ml, catering to all your specific storage needs.

Environmental Responsibility

We're committed to protecting the planet. Our unique and innovative packaging is 100% recyclable and eco-friendly, without compromising on quality or durability.

If you require any more information on this product, please contact our customer service team.