

# PRODUCT DATA SHEET



## **CV0008**

Specimen Tubes, With Unattached Push In, Polythene Stoppers



*Product comes with unattached Push In, Polythene Stoppers*

## **Specifications**

<b>Vial Material</b>	Glass, Soda <i>(Not recommended for autoclave)</i>	●
<b>Outer Diameter</b>	25mm	
<b>Height</b>	50mm	
<b>Closure Type</b>	Push In, Polythene Stoppers	
<b>Closure Liner</b>	Unlined	
<b>Pack Quantity</b>	100	
<b>Packaging</b>	Cardboard Box & Plastic Bag	
<b>Pack Length</b>	400mm	
<b>Pack Width</b>	235mm	
<b>Pack Height</b>	54mm	
<b>Pack Weight</b>	1.3Kg	
<b>Case Quantity</b>	10packs (1000vials)	
<b>Case Length</b>	490mm	
<b>Case Width</b>	420mm	
<b>Case Height</b>	300mm	
<b>Case Weight</b>	13.4Kg	
<b>Pallet Quantity</b>	300packs (30000Vials)	

## **What's in the box?**



Streamline your workflow and simplify your storage with our pack of 100 vials and caps. They neatly stack and nest together, maximizing your space and allowing convenient access whenever needed.

### **Unrivalled Quality**

Each vial is crafted using premium-quality, chemically-resistant glass that ensures excellent durability and protection, and providing unparalleled clarity.

### **Range of Sizes**

We understand that every project is unique. That's why our new range is available in a variety of sizes, from 10mm to 100mm, catering to all your specific storage needs.

### **Environmental Responsibility**

We're committed to protecting the planet. Our unique and innovative packaging is 100% recyclable and eco-friendly, without compromising on quality or durability.

*If you require any more information on this product, please contact our customer service team.*