

# PRODUCT DATA SHEET



## **CV0009**

Specimen Tubes, With Unattached Push In, Polythene Stoppers



*Product comes with unattached Push In, Polythene Stoppers*

## **Specifications**

<b>Vial Material</b>	Glass, Soda ( <i>Not recommended for autoclave</i> )	●
<b>Outer Diameter</b>	19mm	
<b>Height</b>	50mm	
<b>Closure Type</b>	Push In, Polythene Stoppers	
<b>Closure Liner</b>	Unlined	
<b>Pack Quantity</b>	100	
<b>Packaging</b>	Cardboard Box & Plastic Bag	
<b>Pack Length</b>	335mm	
<b>Pack Width</b>	185mm	
<b>Pack Height</b>	55mm	
<b>Pack Weight</b>	1Kg	
<b>Case Quantity</b>	10packs (1000vials)	
<b>Case Length</b>	400mm	
<b>Case Width</b>	360mm	
<b>Case Height</b>	300mm	
<b>Case Weight</b>	9.6Kg	
<b>Pallet Quantity</b>	288packs (28800Vials)	

## **What's in the box?**



Streamline your workflow and simplify your storage with our pack of 100 vials and caps. They neatly stack and nest together, maximizing your space and allowing convenient access whenever needed.

## **Unrivalled Quality**

Each vial is crafted using premium-quality, chemically-resistant glass that ensures excellent durability and protection, and providing unparalleled clarity.

## **Range of Sizes**

We understand that every project is unique. That's why our new range is available in a variety of sizes, from 10mm to 100mm, catering to all your specific storage needs.

## **Environmental Responsibility**

We're committed to protecting the planet. Our unique and innovative packaging is 100% recyclable and eco-friendly, without compromising on quality or durability.

*If you require any more information on this product, please contact our customer service team.*