

# PRODUCT DATA SHEET



## **CV0016**

Screw Neck Vial, Tall Form, With Unattached PP Closure



*Product comes with unattached PP Screw Cap*

## **Specifications**

|                        |                                     |
|------------------------|-------------------------------------|
| <b>Vial Material</b>   | Glass Borosilicate 3 (Autoclavable) |
| <b>Capacity</b>        | 3.5ml                               |
| <b>Outer Diameter</b>  | 12.5mm                              |
| <b>Height</b>          | 46mm                                |
| <b>Closure Type</b>    | Screw, Black Polypropylene          |
| <b>Closure Liner</b>   | Aluminium                           |
| <b>Pack Quantity</b>   | 100                                 |
| <b>Packaging</b>       | Cardboard Box & Plastic Bag         |
| <b>Pack Length</b>     | 210mm                               |
| <b>Pack Width</b>      | 130mm                               |
| <b>Pack Height</b>     | 55mm                                |
| <b>Pack Weight</b>     | 0.4Kg                               |
| <b>Case Quantity</b>   | 16packs (1600vials)                 |
| <b>Case Length</b>     | 440mm                               |
| <b>Case Width</b>      | 280mm                               |
| <b>Case Height</b>     | 240mm                               |
| <b>Case Weight</b>     | 6.5Kg                               |
| <b>Pallet Quantity</b> | 600packs (60000Vials)               |

## **What's in the box?**



Streamline your workflow and simplify your storage with our pack of 100 vials and caps. They neatly stack and nest together, maximizing your space and allowing convenient access whenever needed.

## **Unrivalled Quality**

Each vial is crafted using premium-quality, chemically-resistant glass that ensures excellent durability and protection, and providing unparalleled clarity.

## **Range of Sizes**

We understand that every project is unique. That's why our new range is available in a variety of sizes, from 2 ml to 100 ml, catering to all your specific storage needs.

## **Environmental Responsibility**

We're committed to protecting the planet. Our unique and innovative packaging is 100% recyclable and eco-friendly, without compromising on quality or durability.

*If you require any more information on this product, please contact our customer service team.*