

# PRODUCT DATA SHEET



## **CV0020**

Screw Neck Vial, Squat Form, With Unattached PP Closure



*Product comes with unattached PP Screw Cap*

## **Specifications**

<b>Vial Material</b>	Glass Borosilicate 3 (Autoclavable)
<b>Capacity</b>	28.25ml
<b>Outer Diameter</b>	27.5mm
<b>Height</b>	72mm
<b>Closure Type</b>	Screw, Black Polypropylene
<b>Closure Liner</b>	Aluminium
<b>Pack Quantity</b>	100
<b>Packaging</b>	Cardboard Box & Plastic Bag
<b>Pack Length</b>	445mm
<b>Pack Width</b>	265mm
<b>Pack Height</b>	80mm
<b>Pack Weight</b>	2.3Kg
<b>Case Quantity</b>	3packs (300vials)
<b>Case Length</b>	460mm
<b>Case Width</b>	280mm
<b>Case Height</b>	260mm
<b>Case Weight</b>	7.5Kg
<b>Pallet Quantity</b>	126packs (12600vials)

## **What's in the box?**



Streamline your workflow and simplify your storage with our pack of 100 vials and caps. They neatly stack and nest together, maximizing your space and allowing convenient access whenever needed.

## **Unrivalled Quality**

Each vial is crafted using premium-quality, chemically-resistant glass that ensures excellent durability and protection, and providing unparalleled clarity.

## **Range of Sizes**

We understand that every project is unique. That's why our new range is available in a variety of sizes, from 2 ml to 100 ml, catering to all your specific storage needs.

## **Environmental Responsibility**

We're committed to protecting the planet. Our unique and innovative packaging is 100% recyclable and eco-friendly, without compromising on quality or durability.

*If you require any more information on this product, please contact our customer service team.*