PRODUCT DATA SHEET



<u>CV0021</u>

Screw Neck Vial, Tall Form, With Unattached PP Closure



Product comes with unattached PP Screw Cap

Specifications

Vial Material	Glass Borosilicate 3 (Autoclavable)
Capacity	21.25ml
Outer Diameter	21mm
Height	86mm
Closure Type	Screw, Black Polypropylene
Closure Liner	Aluminium
Pack Quantity	100
Packaging	Cardboard Box & Plastic Bag
Pack Length	452mm
Pack Width	260mm
Pack Height	88mm
Pack Weight	2.7Kg
Case Quantity	6packs (600vials)
Case Length	460mm
Case Width	350mm
Case Height	310mm
Case Weight	15.9Kg
Pallet Quantity	100packs (10000Vials)

What's in the box?



Streamline your workflow and simplify your storage with our pack of 100 vials and caps. They neatly stack and nest together, maximizing your space and allowing convenient access whenever needed.

Unrivalled Quality

Each vial is crafted using premium-quality, chemically-resistant glass that ensures excellent durability and protection, and providing unparalleled clarity.

Range of Sizes

We understand that every project is unique. That's why our new range is available in a variety of sizes, from 2 ml to 100 ml, catering to all your specific storage needs.

Environmental Responsibility

We're committed to protecting the planet. Our unique and innovative packaging is 100% recyclable and eco-friendly, without compromising on quality or durability.