# **PRODUCT DATA SHEET**



# <u>CV0028</u>

Universal Vial, Moulded, Aluminium Closure



## Product comes with fitted aluminium closure

# **Specifications**

Vial Material	Soda (Not recommended for autoclave)
Capacity	30ml
Outer Diameter	28mm
Height	83mm
Closure Type	Screw, Aluminium
Closure Liner	Silica (Please enquire for additional options)
Pack Quantity	144vials
Packaging	Cardboard Box
Pack Length	293mm
Pack Width	265mm
Pack Height	205mm
Pack Weight	7.7Kg
Case Quantity	2packs (288vials)
Case Length	550mm
Case Width	310mm
Case Height	220mm
Case Weight	15.5Kg
Pallet Quantity	68packs (9792Vials)

# What's in the box?



The Universal Vial is a high-quality, moulded vial designed for the secure storage and transport of liquids and small samples. Each vial is equipped with an aluminium closure to ensure a tight seal, maintaining the integrity and sterility of the contents.

## **Available Sizes**

- **7ml Vial CV0029:** Ideal for smaller sample volumes, the 7ml Bijou Vial is compact and convenient for a variety of laboratory applications.
- **30ml Vial CV0028:** The 30ml Universal Vial offers a larger capacity, suitable for storing larger quantities of samples or reagents.

## Key Features

- **Moulded Construction:** The vials are moulded to precise specifications, ensuring uniformity and high performance in demanding laboratory environments.
- **Aluminium Closure:** Each vial is fitted with an aluminium closure, providing a secure, leak-proof seal that protects the contents from contamination and evaporation.
- **Durability:** Made from high-quality materials, these vials are resistant to breakage and are designed to withstand rigorous handling and storage conditions.
- **Versatility:** Suitable for a wide range of applications including clinical, pharmaceutical, and research laboratories.

## **Environmental Responsibility**

We're committed to protecting the planet. Our unique and innovative packaging is 100% recyclable and eco-friendly, without compromising on quality or durability.