# **PRODUCT DATA SHEET**



### CV0032

Rolled Rim Vials, With Unattached Clip On / Snap On Closure



Product comes with unattached PE clip on / snap on closure

## **Specifications**

Vial Material	Glass Borosilicate 3 (Autoclavable)
Capacity	28.25ml
Outer Diameter	27.5mm
Height	72mm
Closure Material	Polyethylene (PE)
Closure Type	Clip On / Snap
Closure Liner	N/A
Pack Quantity	100
Packaging	Cardboard Box & Plastic Bag
Pack Length	450mm
Pack Width	260mm
Pack Height	75mm
Pack Weight	1.9Kg
Case Quantity	5packs (500vials)
Case Length	470mm
Case Width	280mm
Case Height	410mm
Case Weight	10.5Kg
Pallet Quantity	200packs (20000vials)

#### What's in the box?



Streamline your workflow and simplify your storage with our pack of 100 vials and caps. They neatly stack and nest together, maximizing your space and allowing convenient access whenever needed.

#### **Unrivalled Quality**

Each vial is crafted using premium-quality, chemically-resistant glass that ensures excellent durability and protection, and providing unparalleled clarity.

#### **Range of Sizes**

We understand that every project is unique. That's why our new range is available in a variety of sizes, from 7ml to 30ml, catering to all your specific storage needs.

#### **Environmental Responsibility**

We're committed to protecting the planet. Our unique and innovative packaging is 100% recyclable and eco-friendly, without compromising on quality or durability.

If you require any more information on this product, please contact our customer service team.