

PRODUCT DATA SHEET



CV0038

Specimen Tubes, With Unattached Push In, Polythene Stoppers



Product comes with unattached Push In, Polythene Stoppers

Specifications

Vial Material	Glass, Soda <i>(Not recommended for autoclave)</i>	●
Outer Diameter	25mm	
Height	100mm	
Closure Type	Push In, Polythene Stoppers	
Closure Liner	Unlined	
Pack Quantity	100	
Packaging	Cardboard Box & Plastic Bag	
Pack Length	400mm	
Pack Width	235mm	
Pack Height	104mm	
Pack Weight	2.7Kg	
Case Quantity	10packs (1000vials)	
Case Length	470mm	
Case Width	280mm	
Case Height	410mm	
Case Weight	28Kg	
Pallet Quantity	120packs (12000Vials)	

What's in the box?



Streamline your workflow and simplify your storage with our pack of 100 vials and caps. They neatly stack and nest together, maximizing your space and allowing convenient access whenever needed.

Unrivalled Quality

Each vial is crafted using premium-quality, chemically-resistant glass that ensures excellent durability and protection, and providing unparalleled clarity.

Range of Sizes

We understand that every project is unique. That's why our new range is available in a variety of sizes, from 10mm to 100mm, catering to all your specific storage needs.

Environmental Responsibility

We're committed to protecting the planet. Our unique and innovative packaging is 100% recyclable and eco-friendly, without compromising on quality or durability.

If you require any more information on this product, please contact our customer service team.