

PRODUCT DATA SHEET



CV0044

Scintillation Vial, Unattached White Screw Cap



Product comes with unattached PP Screw Cap

Specifications

Vial Material	Glass Borosilicate 3 (Autoclavable)
Capacity	20ml
Outer Diameter	27mm
Height	51mm
Closure Type	Screw, White Polypropylene
Closure Liner	Foil Faced Liner
Pack Quantity	100
Packaging	Cardboard Box & Plastic Bag
Pack Length	447mm
Pack Width	261mm
Pack Height	61mm
Pack Weight	1.9Kg
Case Quantity	5packs (500vials)
Case Length	470mm
Case Width	280mm
Case Height	340mm
Case Weight	9.4Kg
Pallet Quantity	200packs (20000Vials)

What's in the box?



Streamline your workflow and simplify your storage with our pack of 100 vials and caps. They neatly stack and nest together, maximizing your space and allowing convenient access whenever needed.

Unrivalled Quality

Each vial is crafted using premium-quality, chemically-resistant glass that ensures excellent durability and protection, and providing unparalleled clarity.

Environmental Responsibility

We're committed to protecting the planet. Our unique and innovative packaging is 100% recyclable and eco-friendly, without compromising on quality or durability.

If you require any more information on this product, please contact our customer service team.