

PRODUCT DATA SHEET



CV0045

Bijou, Tubular, White Polypropylene Closure




Product comes with unattached polypropylene closure

Specifications

Vial Material	Glass, Soda <i>(Not recommended for autoclave)</i>	●
Capacity	7ml	
Outer Diameter	20mm	
Height	43mm	
Closure Type	Screw, White Polypropylene	
Closure Liner	Unlined	
Pack Quantity	100vials	
Packaging	Cardboard Box	
Pack Length	408mm	
Pack Width	200mm	
Pack Height	49mm	
Pack Weight	1Kg	
Case Quantity	10packs (1000vials)	
Case Length	500mm	
Case Width	230mm	
Case Height	290mm	
Case Weight	10.5Kg	
Pallet Quantity	270packs (27000vials)	

What's in the box?



The Bijou Vial, tubular vial designed for the secure storage and transport of liquids and small samples. Each vial is equipped with a polypropylene closure, maintaining the integrity and sterility of the contents.

Available Sizes

- **7ml Vial CV0045:** Ideal for smaller sample volumes, the 7ml Bijou Vial is compact and convenient for a variety of laboratory applications.
- **28ml Vial CV0046:** The 28ml Universal Vial offers a larger capacity, suitable for storing larger quantities of samples or reagents.

Key Features

- **Tubular Construction:** The vials are produced to precise specifications, ensuring uniformity and high performance in demanding laboratory environments.
- **Polypropylene Closure:** Each vial has an unattached polypropylene closure, providing a secure seal that protects the contents from contamination and evaporation.
- **Versatility:** Suitable for a wide range of applications including clinical, pharmaceutical, and research laboratories.

Environmental Responsibility

We're committed to protecting the planet. Our unique and innovative packaging is 100% recyclable and eco-friendly, without compromising on quality or durability.

If you require any more information on this product, please contact our customer service team.